

VELOCITY
orthopedics

**Value
Through
Innovation**

CleanPass™

Velocity®...

A premium line of
instrumentation for
arthroscopic surgery.

Ideal for arthroscopic labral
and rotator cuff repairs, the
CleanPass™ Suture Passer
integrates best-in-class
features into an elegant, low
profile device.

*Tip configurations to
suit the most demanding
applications*

**Versatile, Robust,
Cost-Effective**

Versatile, Robust, Cost-Effective


CleanPass™

- Secure-Lok™ suture needle connection - robust and reliable
- Reusable handle and semi-reusable suture hooks - consistent performance at a fraction of the disposable cost
- Two convenient suture loading ports for ergonomic flexibility
- Bevel ground, sharp, 1.6mm needles - atraumatic suture passage steadied by a robust hub
- Fits easily through an 8mm cannula
- Tip configurations to suit the most demanding applications
- Monofilament suture shuttle available
- Roller wheel replacement kit available for easy field replacement

Secure-Lok. Line up notch, tighten securely



Two Loading Ports, top and rear



SUTURE PASSER ORDERING INFORMATION		
4100	CleanPass™ Suture Passing Handpiece	
4101	CleanPass™ Suture Passing Needle - Crescent	
4102	CleanPass™ Suture Passing Needle - 45 Left	
4103	CleanPass™ Suture Passing Needle - 45 Right	
4104	CleanPass™ Suture Passing Needle - 90 Up	
4105	CleanPass™ SI Roller Replacement Kit	
4106	CleanPass™ Suture Shuttle	

4050 TRAY ORDERING INFORMATION		
4100	CleanPass™ Suture Passing Handpiece	
4101	CleanPass™ Suture Passing Needle - Crescent	
4102	CleanPass™ Suture Passing Needle - 45 Left	
4103	CleanPass™ Suture Passing Needle - 45 Right	
4104	CleanPass™ Suture Passing Needle - 90 Up	
4105	CleanPass™ SI Roller Replacement Kit	
4106	CleanPass™ Suture Shuttle	
4032	Suture Cutter	
4006	Suture Grasper	
4004	Retriever	
4008	Javelin - 30	



Velocity® Orthopedics, Inc.
 10700 Jersey Blvd., Suite 360
 Rancho Cucamonga, CA 91730 USA
 Tel: 909.987.4343 Fax: 909.987.5800
 © 2017 Velocity Orthopedics, Inc.

velocityorthopedics.com

Made in USA

P/N MM006 2/17