



## VP150 Vista Arthroscopy Pump for arthroscopy

- Minimized dimensions
- Optimized performance
- Automatic instrument detection

# VP150 Vista Arthroscopy

## Pump for arthroscopy

- 1.5 l/min max. flow
- Optimized flow regulation done by automatic instrument detection and software control
- Wash function to improve visibility within the joint
- Flexible use due to pole mounting
- Highest patient safety due to transponder protected tube set

### Device

Item No.	Description
VP150	Pump for arthroscopy

### Accessories

Part No.	Description
	Tube set with Spike and Luer-Lock connector male:
VP150-10	Disposable, only Inflow, transponder protection (1 packaging = 10 pcs.)



Pole mounting

## Tube sets for best performance and patient comfort

- Safe and reliable tube set connection
  - Ensures the optimal performance
- RFID Transponder Technology
  - Minimizes the potential risk of misuse
  - High protection against copies



Tube set with Spike and Luer-Lock connector male, disposable, only Inflow, transponder protection

#### Technical Data

Pressure range	5-150 mmHg
Max. flow	1.5 l/min $\pm$ 10 %
Mains	Medical world power supply, 100-240 V AC, 50-60 Hz
Power consumption	60 VA
Current consumption	100 V: 1.2 A 240 V: 0.6 A
Fuse designation	2x T 3.15 AH, 250 V
Classifications	Protection class I, Type BF, IP 41
Dimensions (w x h x d)*	240 x 103 x 221 mm
Weight	2.8 kg

\*without toggle screw

Distributor:

**OrthoTec International**

7440 SW 50 Terrace, # 108

Miami, Fl. 33155 – USA

Phone +1-305-662-2855

All technical data in our other sales information are subject to modifications. Any liability that the data and illustrations are complete, actual or correct is excluded.

Products may not available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact us if you have questions about the availability of our products in your area.

All images, text and material are protected by copyright and may not be copied, printed or reproduced without written consent from the copyright owner.